


# nAspro-bby

Blue Blocking Yellow

## Nano Aspheric Precision Optics

Aspheric Foldable Posterior Chamber  
Intraocular lens



Enhanced 360° Square Edge

Blue Blocking Optics

CE 2460

 **APPASAMY  
ASSOCIATES**  
Empowering Vision®

# nAspro-bby

Blue Blocking Yellow

## Technical Specification

Model	NASY 207
Overall length	12.50 mm
Optic Diameter	6 mm
Haptic Angulations	0°
Optics Design	Biconvex
Diopter Range	6 D - 30.0 D Diopter (0.5 Diopter)
Material	Hydrophilic Acrylic (yellow color)
Refractive Index	1.461
Delivery System	Disposable Delivery System
A constant	118.0
Special Features	Aspheric Optics with 360° Square edge
Water Content	26%

Manufactured by:

 **APPASAMY OCULAR DEVICES (P) LTD**

R.S.NO.9/1,2&3, N.H.45-A, Villupuram Main Road,  
Vadamangalam, Pondicherry - 605 105, INDIA

Marketed by :

 **IMPLANTEC**

**IMPLANTEC S.A.**

Perdriel 1624 · Ciudad Autónoma de Buenos Aires  
República Argentina

Tel: (+54 11) 43 02 71 11 · [info@implantec.com.ar](mailto:info@implantec.com.ar)

[www.implantec.com.ar](http://www.implantec.com.ar)